



ziehm imaging

ARTICULATING MONITOR ARM

Maximum
Flexibility and
Ergonomic Comfort



- Horizontal monitor rotation: +/- 40°
- Vertical monitor swiveling: +/- 10°
- Horizontal movement range: +/- 115°
- Vertical adjustment range: + 9 in / - 6 in
- Maximum range: 32 in

Position the Remote Vision Center (RVC) and the Position Control Center (PCC) on the back of the monitor cart for easy access. The Articulating Monitor Arm* is available for all Ziehm Vision systems.

* Articulating Monitor Arm is an option for purchase that can be utilized with either the 32" monitor or the 19" monitor

Ziehm Imaging
6280 Hazeltine National Dr | Orlando, FL 32822 | USA
Phone 800 503 4952 | Fax 321 445 5514
mail@ziehm.com | www.ziehm.com



© 2023 Ziehm Imaging, a division of Ziehm-Orthoscan, Inc. All Rights Reserved. Ziehm Imaging is constantly improving its products and reserves the right to change these specifications without notice. 113-0061 Rev. D 07/2023

Nationwide toll-free
**SERVICE & TECHNICAL
PHONE SUPPORT**
866.949.4346

