



ziehm imaging

CONTROL MODULES

Precise and Timesaving



Optional
Mobile Stand

Easily operate and move the motorized C-arm in all four axes and keep track of orientation and object positioning entirely in the sterile field; minimizing the time needed and maximizing precision.

- The **Remote Vision Center** is a touchscreen control panel that can be mounted directly to an OR table or on an optional mobile stand and provides the same functions as the touchscreen control panels on the monitor cart and C-arm
- The **Position Control Center (joysticks)** is used by the operator to move the C-arm into the exact desired position and with the freely selectable isocenter, any given anatomical structure can be displayed from different angles without having to re-adjust the C-arm
- Easily save and recall up to **3 positions**, for example to switch between AP, lateral and oblique positions

The Remote Vision Center and the Position Control Center are optional control panels that can be connected to the monitor cart or used on the optional mobile stand for the Ziehm Vision C-arms.

Ziehm Imaging
6280 Hazeltine National Dr | Orlando, FL 32822 | USA
Phone 800 503 4952 | Fax 321 445 5514
mail@ziehm.com | www.ziehm.com



© 2023 Ziehm Imaging, a division of Ziehm-Orthoscan, Inc. All Rights Reserved. Ziehm Imaging is constantly improving its products and reserves the right to change these specifications without notice. 113-0009 Rev. E 07/2023

Nationwide toll-free
**SERVICE & TECHNICAL
PHONE SUPPORT**
866.949.4346

